


# SMART INSULIN PENS: NovoPen® 6 and NovoPen Echo® Plus How to upload data with Glooko®

1. You will need a compatible smartphone with an NFC chip to upload data from a smart insulin pen.

2. Start Glooko® Mobile App and log in. The App is available on [App Store](#) (iOS) and [Google Play](#) (Android).



3. Make sure that the pen's dose button is fully pushed, and the dose counter is at zero.

4.

(new Glooko user)

(existing Glooko user)

(first time connecting device)

(after device has been connected)

On the Home tab, click 'Sync' in the top right corner.

If connecting your device for the first time, click 'New Device' in the top right and select 'Smart Pens'. If you have already connected your device, select 'Sync' next to the device you wish to upload.

5. Follow on-screen instructions to complete the upload.

6. After successfully uploading, you can view your insulin data in a variety of charts and reports.

## HOW TO LOCATE THE NFC CHIP ON YOUR PHONE

Normally the NFC chip is fitted at the centre of the back of the phone, but this can vary between different models. Please contact the mobile phone manufacturer if there are any connectivity issues between the pen and the phone and to find out where the NFC chip is located on your phone.

NovoPen® 6 and NovoPen Echo® Plus are registered trademarks owned by Novo Nordisk.